

## T3976 Tri-Cornered Beakers

### I. APPLICATIONS

1. Transfer of chemicals and other liquids

### II. FEATURES

1. Materials: Strong, medical-grade polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. Transparent wall permits easy viewing of tube contents
4. Autoclavable

### III. BEAKER DESIGN

1. Unique three-corner design
2. Available in 50ml,100ml, 250ml, 400ml, 800ml and,1000ml
3. Molded graduations
4. Clearly visible with either clear or turbid samples

### IV. PACKAGING

1. Bulk Packaging, 25 pieces/sleeve and 4 sleeves/case

### V.DESCRPTION

Item	Size	Packaging	Sterile
93-0050	50ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0100	100ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0250	250ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0400	400ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-0800	800ml Beaker	25 pieces/sleeve, 4 sleeves/case	No
93-1000	1,000ml	25 pieces/sleeve, 4 sleeves/case	No

### VI. STERILITY

1. N/A

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

### VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

### IX. PRODUCT PARAMETER

Dimensions=mm

Cat # <sup>⊖</sup>	Description <sup>⊖</sup>	A <sup>⊖</sup>	B <sup>⊖</sup>	C <sup>⊖</sup>	D <sup>⊖</sup>	E <sup>⊖</sup>	F <sup>⊖</sup>	G <sup>⊖</sup>
93-0050	Beaker <sup>⊖</sup>	67 <sup>⊖</sup>	58 <sup>⊖</sup>	O.D.	O.D. 43.1 <sup>⊖</sup>	O.D. 46.9 <sup>⊖</sup>	31 <sup>⊖</sup>	57 <sup>⊖</sup>
93-0100	Beaker <sup>⊖</sup>	84 <sup>⊖</sup>	72 <sup>⊖</sup>	O.D.	O.D. 53.7 <sup>⊖</sup>	O.D. 58.4 <sup>⊖</sup>	39 <sup>⊖</sup>	71 <sup>⊖</sup>
93-0250	Beaker <sup>⊖</sup>	101 <sup>⊖</sup>	88 <sup>⊖</sup>	O.D.	O.D. 68 <sup>⊖</sup>	O.D. 74 <sup>⊖</sup>	49 <sup>⊖</sup>	88 <sup>⊖</sup>
93-0400	Beaker <sup>⊖</sup>	114.7	99.9 <sup>⊖</sup>	O.D. 52 <sup>⊖</sup>	O.D. 82 <sup>⊖</sup>	O.D. 84 <sup>⊖</sup>	55.6 <sup>⊖</sup>	110.3
93-0800	Beaker <sup>⊖</sup>	133.1	115.2	O.D. 65 <sup>⊖</sup>	O.D. 104 <sup>⊖</sup>	O.D. 105.1 <sup>⊖</sup>	65 <sup>⊖</sup>	137 <sup>⊖</sup>
93-1000	Beaker <sup>⊖</sup>	140.6	121.7	O.D.	O.D. 112.7 <sup>⊖</sup>	O.D. 114.1 <sup>⊖</sup>	70.9 <sup>⊖</sup>	143 <sup>⊖</sup>

